
VISULAS YAG III

Gebrauchsanweisung
User's Manual
Mode d'emploi
Instrucciones para el uso

ZEISS

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Knowledge of this manual is required for operation of the instrument. You should therefore familiarise yourself with the contents of this manual and pay special attention to instructions concerning the safe operation of the instrument.

We reserve the right to make changes to the product in light of technical developments; this user manual will not be updated.

© Unless expressly authorised, dissemination or duplication of this document and commercial exploitation or communication of its contents is not permitted. Persons in contravention of this copyright are liable to pay compensation for damages.
All rights reserved in the events of granting of patents or registration as a utility patent.

	Page
Copyright.....	1
Contents / Figures.....	3
Safety Notes.....	6
General	6
Standards and regulations	7
Installation and operating instructions	8
Laser safety	11
Instrument safety systems.....	15
Warning and information labels.....	16
Description	18
Design of the VISULAS YAG III.....	18
Beam delivery in the YAG laser head	19
Design of the VISULAS YAG III Combi.....	20
Beam delivery in the LSL YAG III Combi laser slit lamp.....	21
Controls	22
VISULAS YAG III control unit	22
VISULAS 532s control panel for VISULAS YAG III Combi .	24
Laser slit lamp	26
Instrument base with headrest	28
Intended use	30
Contact lenses.....	30
Principle of Nd:YAG laser treatment	31
Contraindications	32
Operation	33
Installation and transport of the VISULAS YAG III.....	33
Installation of the VISULAS YAG III Combi.....	35
Moving the IT 4L instrument table with the instrument...	37
Laser treatment using the VISULAS YAG III /	
VISULAS YAG III Combi	38
Adjusting slit lamp and eyepieces	38
Start-up and laser treatment	38
VISULAS YAG III: Description of menus.....	39
Menu navigation.....	39
Power-up mode	40
Diagnosis mode	41
Treatment STANDBY mode	42
Information mode	44
Setting the system time	46
Treatment READY mode	47
VISULAS YAG III Combi: Description of menus.....	48
Menu navigation.....	48
VISULAS 532s power-up mode.....	49
Diagnosis mode	50
VISULAS YAG III power-up mode	51
YAG III Treatment mode.....	51
STANDBY	52
READY	54
Using the aiming beam system with the VISULAS YAG III /	
VISULAS YAG III Combi	55
Focus shift and defocusing	56

Maintenance / Miscellaneous.....	57
Troubleshooting tables.....	57
Error messages on the VISULAS YAG III control unit	57
Faults during operation of the VISULAS YAG III Combi.....	58
Error messages and information on the VISULAS 532s control unit.....	58
Laser slit lamp faults	60
Replacing the fuse on the VISULAS YAG III control unit	61
Replacing the halogen lamp in the laser slit lamp	62
Care of the instrument.....	63
Cleaning optical components.....	63
Cleaning of painted surfaces.....	63
Disposal	64
Disposal of the product within the EU	64
Safety checks	65
Calibration of energy measurement	66
Requirements	66
Calibration procedure	67
Ordering data, accessories and spare parts	69
Technical Specifications	70
VISULAS YAG III laser system.....	70
LSL YAG III laser slit lamp	71
LSL YAG III Combi laser slit lamp	71
Ambient conditions for overall system	72
Manufacturer's Declaration.....	73
Fig. 1	Warning and information labels..... 16
Fig. 2	Location of warning and information labels on the VISULAS YAG III..... 17
Fig. 3	Location of warning and information labels on the LSL YAG III Combi laser slit lamp..... 17
Fig. 4	VISULAS YAG III instrument components
Fig. 5	LSL YAG III beam delivery..... 19
Fig. 6	VISULAS YAG III Combi instrument components
Fig. 7	Beam delivery in the LSL YAG III Combi laser slit lamp..... 21
Fig. 8	VISULAS YAG III control unit
Fig. 9	Control unit/connectors
Fig. 10	Removing the VISULAS 532s control panel..... 25
Fig. 11	VISULAS YAG III laser slit lamp
Fig. 12	LSL YAG III Combi laser slit lamp
Fig. 13	Instrument base with headrest..... 29
Fig. 14	Principle of photodisruption using the Nd:YAG laser
Fig. 15	VISULAS YAG III in transport case
Fig. 16	IT 4L with table mounting..... 35
Fig. 17	IT 4L with laser consoles
Fig. 18	Binocular tube and eyepieces..... 36
Fig. 19	VISULAS YAG III menu navigation
Fig. 20	Power-up mode..... 40
Fig. 21	Diagnosis mode
Fig. 22	Treatment STANDBY mode..... 42
Fig. 23	Information mode..... 44

Fig. 24	Setting the system time	46
Fig. 25	Treatment READY mode	47
Fig. 26	VISULAS YAG III Combi menu navigation	48
Fig. 27	Power-up mode	49
Fig. 28	Diagnosis mode	50
Fig. 29	System test	51
Fig. 30	Combination mode activated	51
Fig. 31	Main menu after selecting the YAG laser as applicator ...	51
Fig. 32	Treatment STANDBY mode	52
Fig. 33	Treatment READY mode	54
Fig. 34	Focusing the aiming beam	55
Fig. 35	Distorted aiming beam	55
Fig. 36	Focus shift control	56
Fig. 37	Replacing the fuse	61
Fig. 38	Replacing the halogen lamp	62

General

This instrument has been developed and tested in accordance with national and international regulations. A high degree of instrument safety is thus ensured.

This chapter contains a summary of the most important safety-relevant information.

All safety notes and information in this manual and on the labels on the instrument must be observed. These notes are marked by pictograms.

	Warning	Danger to the user!
	Caution	Danger to the instrument!
	Symbol on type label	Type B medical device in accordance with IEC 60601-1
		Disconnect the power cable before opening the device!
		Caution: Laser radiation!
	Note	Information and notes for a better understanding of the user manual

**Note for VISULAS YAG III Combi:**
Special note when using the VISULAS YAG III Combi.

**Warning**
Correct operation of the device is imperative for its safe function. You should therefore ensure that you are thoroughly familiar with this user manual before using the instrument.

Laser damage may occur if the instrument is operated or used for the treatment of patients in a manner other than that described in this user manual.

Operating instructions for other equipment used with this instrument must be observed.

Further information may be obtained from Carl Zeiss Meditec authorised specialists or from authorised representatives.

Standards and regulations

- ❑ Carl Zeiss Meditec works according to a certified quality management system.
- ❑ In accordance with the relevant standards, this instrument is equipped with an energy display, a key switch, an interlock device and all required warning and information signs.
- ❑ All relevant national accident prevention regulations for this instrument must be observed.
- ❑ In some countries, national regulations specify that this instrument may only be used under the supervision of a medical doctor.

Warning

In the USA, this device may only be sold to a medical doctor or on the basis of an order from a medical doctor.

- ❑ This instrument is a class IV (4) laser device. The safety regulations applying to this laser classification must be observed
- ❑ This instrument is a class IIb medical instrument as defined by the European Medical Device Directive (MDD).
- ❑ The instrument complies with the EU Medical Device Directive 93/42/EEC and its national implementation in the form of the German Medical Products Act (MPG) (→ *Manufacturer's Declaration* on page 73).
- ❑ To ensure proper function of the instrument, regular safety inspections should be carried out. The instrument should be inspected by a Carl Zeiss Meditec authorised specialist annually and the results entered in the instrument logbook. See also chapter *Safety checks* of this user manual.
- ❑ An instrument logbook must be kept.

Installation and operating instructions

General

- ❑ This instrument may only be operated by appropriately qualified and trained personnel. It is the responsibility of the institution operating the instrument to train and instruct operating personnel.
- ❑ People working in the laser hazard area must be given instruction concerning safety regulations, safety precautions and operation of the instrument at least once a year. A written record of this instruction, listing the participants, must be kept.
- ❑ The user manual and instrument logbook should be at hand for persons operating the instrument at all times.
- ❑ The instrument should be used for the applications described only. Damage occurring through the instrument being used in a manner other than that intended is the responsibility of the user.
- ❑ Do not operate the instruments supplied
 - in potentially explosive atmospheres,
 - in the presence of volatile anaesthetics or flammable solvents, such as alcohol, benzene, etc.
- ❑ Do not store or use this device in damp rooms. Do not expose the instrument to water splashes, dripping water or sprayed water.
- ❑ Do not place fluid-filled containers on top of the control unit.
- ❑ Switch the instrument off immediately if smoke, sparks or unusual noises are observed. Do not use the instrument until it has been repaired by a Carl Zeiss Meditec authorised specialist.
- ❑ If, after switching on the instrument with the key switch, the display fails to light up, shut down the instrument, disconnect the power cable and attach a sign which clearly states that the instrument is out of order. Do not use the instrument until it has been repaired by a Carl Zeiss Meditec authorised specialist.
- ❑ Do not use force when inserting connectors. If the connector cannot readily be inserted into the jack or socket, make sure that you are using the correct connector for the socket. Do not pull the cables. Defective connectors should be repaired by a Carl Zeiss Meditec authorised specialist.
- ❑ Modifications and repairs to this instrument and any devices operated in combination with the VISULAS 532s are to be carried out by Carl Zeiss Meditec authorised specialists only. The manufacturer accepts no liability for damages caused by unauthorised tampering with the instrument. Such actions will also render any warranty claims invalid.
- ❑ This instrument may only be used with accessory equipment supplied by Carl Zeiss Meditec.

- ❑ For safety reasons and to ensure adherence to the guarantee conditions, the laser console may only be opened by specialist personnel with express written authorisation from Carl Zeiss Meditec.
- ❑ Disconnect the system from the power supply before opening the instrument, carrying out servicing or changing the fuses.
- ❑ Do not use mobile phones or other devices not compliant with EMC class B in the vicinity of this equipment, as their signals may cause the equipment to malfunction. The effect of radio signals on medical instruments is dependent on various factors and is therefore not predictable. To avoid electromagnetic interference, the instrument must be installed and operated as described in the user manual and using only those components supplied by Carl Zeiss Meditec.
- ❑ Do not connect the system via additional portable multiple sockets or extension cables.
- ❑ The manufacturer accepts no liability for damages arising as a result of unauthorised access to the interior of the instrument or the use of unauthorised accessories.
Such actions will also render any warranty claims invalid.

Warning

Internal components remain electrically live after the instrument has been switched off using the key switch.

To completely disconnect from the power supply, the instrument must be unplugged.

Requirements for operation of the instrument

The instrument will be installed by a Carl Zeiss Meditec authorised specialist. Please ensure that the following requirements for operation are always fulfilled.

- ❑ The ambient conditions for proper use are maintained (see *Technical Specifications* on page 70).
- ❑ The instrument is plugged into a power outlet, which has a properly connected protective earth connection.
- ❑ The power cable intended for use with this instrument is used.
- ❑ No external damage to the instrument is visible.
- ❑ All cables and connectors are in perfect working condition.
- ❑ Pay special attention to warning and information labels and any parts, such as screws and surfaces on the instrument painted red.
- ❑ Do not cover ventilation grids or slits.

When operating the instrument, always

- ☐ Use the lowest possible energy necessary to achieve the desired physiological effect.
- ☐ Before each treatment, check that the focus shift control has been set to the proper position (see *Focus shift and defocusing*, p. 56).
- ☐ Focus the aiming beam as carefully as possible onto the area to be treated. Never fire a laser pulse if you are unable to see the aiming beam in the target area.
- ☐ Always switch off the aiming beam when the instrument is not in use.
- ☐ Because the aiming beam follows the same optical path through the beam delivery system as the therapy beam, it represents a good method for checking that the beam delivery system is undamaged. If the aiming beam spot at the distal end of the beam delivery system is not visible, is weak or appears diffuse, it may indicate damage to or improper function of the beam delivery system.

After use

- ☐ Always use the key switch to switch the instrument off.
- ☐ Always remove the key when the instrument is not in use.
- ☐ Do not allow unauthorised personnel access to the key switch key.

Safe working order

This instrument is a high quality technical product. To ensure optimal performance and that the instrument remains in safe working order, the instrument must be inspected by a Carl Zeiss Meditec authorised specialist at least once a year (see *Safety checks*).

Note for VISULAS YAG III Combi:

Warning

Optical fibres must not be bent sharply, kinked or incorrectly fastened, as this may cause damage to the instrument and/or injury to the patient or user.

Laser safety

Safety precautions for the laser hazard area

This instrument is a class 4 laser device. The institution operating this instrument must take safety precautions to prevent any hazard. Applicable national and international regulations must be observed.

We draw your attention to some important aspects of these regulations:

- ❑ The laser hazard area is the area in which the maximum permissible exposure (MPE) can be exceeded. This includes areas into which the laser beam might unintentionally be deflected by a reflecting surface.
- ❑ The laser hazard area should be as small as possible and should be demarcated using suitable screens. Unauthorised persons should be prevented from entering the laser hazard area. The number of persons in the laser hazard area should be restricted to as few as possible. Where suitable screens are not in use, the entire treatment room must be considered as a laser hazard area.
- ❑ Warning signs must be placed at the entrance to the laser hazard area.
- ❑ A warning lamp must be installed at the entrance indicating when the laser is in operation.
- ❑ The laser hazard area must be clearly demarcated and marked as such during laser operation. A remote interlock contact can be fitted. This switches off the laser automatically if anyone enters the laser hazard area during operation. Further safety precautions should be clarified with local authorities or with your laser safety officer.
- ❑ All objects within the laser hazard area, including the floor, should have diffusely reflecting surfaces or should be covered with diffusely reflecting material.
- ❑ Only the patient and appropriately trained personnel may remain within the laser hazard area. These persons must wear suitable laser safety goggles when the laser is in operation.
- ❑ Persons working in the laser hazard area must be given instruction concerning safety regulations, safety precautions and operation of the instrument at least once a year. A written record of this instruction, listing the participants, must be kept.

Safety officer

The institution operating the instrument and the safety officer are responsible for implementing all safety measures, such that the patient, the attending physician and any persons present during laser treatment are not exposed to any danger. In Germany BGV B2 must be observed. In other countries, the relevant national regulations must be observed.

The responsibilities of the safety officer, who is to be appointed in writing by the institution operating the laser, include:

- ☐ Implementation of safety measures.
- ☐ Instruction of all persons working with the instrument in safety measures and correct operation of the instrument.
- ☐ Marking the laser hazard area.
- ☐ Checking the warning signs and equipment.
- ☐ Correct therapeutic use of the instrument.
- ☐ Safekeeping of the key switch key.
- ☐ Safekeeping of the laser device.
- ☐ Proper reconnection of the device after relocation.
- ☐ Correct keeping of the instrument logbook.

Laser safety goggles

- ☐ All persons present in the laser hazard area during laser operation must wear laser safety goggles. The patient's eye which is not being treated must also be properly protected.
- ☐ Exception:
The physician need not wear safety goggles when observing the area to be treated through the Carl Zeiss Meditec laser slit lamp.
- ☐ The institution operating the instrument is responsible for providing proper laser safety goggles. The safety goggles must be approved for the wavelength of the laser and for the type and intensity of the laser radiation (in Europe, in compliance with DIN EN 207).

VISULAS YAG III: IR 1064 L5

Note for VISULAS YAG III Combi:

Laser safety goggles for operation in combination with the VISULAS 532s: D 450 ... 532 L5

Patient safety

- ❑ The most important parameter in laser treatment is the energy flux density at the treatment area, i.e. the applied laser energy divided by the area of the laser spot.
On the VISULAS YAG III, the energy of the laser pulse can be pre-set.

Note for VISULAS YAG III Combi:

The most important parameter in laser treatment is the power density at the treatment area, i.e. the applied laser power divided by the area of the laser spot. To achieve a constant physiological effect, the power output must also be changed when the spot diameter is changed.

On the VISULAS YAG III Combi, the power output (VISULAS 532s mode) or the energy (VISULAS YAG III mode) can be pre-set.

Explosion and fire hazards

- ❑ Do not use the laser in combination with flammable anaesthetics.
- ❑ Keep explosive materials away from the laser hazard area. Easily flammable materials may cause fire.
- ❑ Laser beams can set fire to many explosive or flammable gases and liquids, including some solvents used in preparing for surgical procedures.
- ❑ Keep flammable drapes, surgical gowns, gauze and other flammable materials out of the beam path. We recommend using non-flammable materials and instruments and wearing gowns and clothing that do not easily catch fire.
- ❑ Keep a fire extinguisher near the laser device.

Hazards caused by direct exposure and reflection

- ❑ The instrument emits invisible and visible laser radiation that is hazardous to the human eye and skin.
- ❑ Do not expose the eye or skin to direct or scattered radiation! Metal and other surfaces reflect laser beams. All persons present in the laser hazard area must therefore wear laser safety goggles as a precaution against unintentional direct or reflected laser radiation.
- ❑ Remove or cover any shiny, reflecting objects within the laser hazard area.
- ❑ Windows and reflecting walls should also be covered with non-flammable cloths.
- ❑ Safety precautions for dealing with noxious gases, dust and vapours, secondary radiation or explosive gas mixtures which may be produced by exposure of materials within the laser hazard area to laser radiation should be put in place.
- ❑ Only medical instruments which largely exclude hazardous reflections due to their design and surface finish should be brought into the laser beam path.

Faults

- ❑ If, during testing or operation of the laser system, faults occur which cannot be remedied by following the instructions given in the section *Maintenance / Miscellaneous*, the instrument should be switched off, the power cable disconnected and the instrument clearly labelled as being out of order.
- ❑ Inform the responsible service technician.

Instrument safety systems

Safety system	Effect
Key switch (1, Fig. 8)	The laser cannot be switched on by unauthorised persons. The user must remove the key after the treatment.
Remote interlock contact (2, Fig. 9)	The ability to fire the laser can be made dependent on the status of an external contact, for example on a door closer. A Carl Zeiss Meditec authorised specialist will be happy to assist, if you wish to connect the laser to a remote interlock. If the instrument is not connected to a remote interlock, a shorting plug which bridges the respective contacts must be inserted in the interlock socket (The device is delivered with the shorting plug inserted.).
Monitoring of therapy beam energy	Laser firing is automatically locked when the metered energy is outside the pre-defined range.
Safety monitoring	The electronic systems monitor various functions and values. If any of these values are outside the permitted range, the system prevents the laser from being fired.
LASER STOP button (4, Fig. 8)	The LASER STOP button is intended for control of unforeseeable hazardous situations. When this button is pressed, the system goes into STANDBY mode and all functions on the control unit are locked. To resume operation, press the LASER STOP button again.
Laser warning lamp (1, Fig. 25)	This lamp lights in READY mode and in STANDBY mode, when the aiming beam is on.

Warning and information labels


**Caution**

All warning and information labels must be observed!


Warning and information labels and their location on the instrument are shown in Fig. 1, Fig. 2 and Fig. 3.

Should it become apparent that one of these labels is missing from your instrument, please contact Carl Zeiss Meditec or a Carl Zeiss Meditec authorised representative.


Position 1


Position 2

Position 3


Position 4


Position 7

Position 5

Position 8

Position 6

Fig. 1 Warning and information labels


Fig. 2 Location of warning and information labels on the VISULAS YAG III



Fig. 3 Location of warning and information labels on the LSL YAG III Combi laser slit lamp

Design of the VISULAS YAG III

The instrument consists of a laser head (1) mounted to a laser slit lamp (2) and a control unit (3).

An optional foot switch (4) is available. The control unit menu is used to switch between the hand and foot switches.

The VISULAS YAG III can be mounted to any table with an even table-top.

- 1 VISULAS YAG III laser head
- 2 LSL YAG III laser slit lamp
- 3 VISULAS YAG III control unit (laser console)
- 4 Foot switch (optional)

Fig. 4 VISULAS YAG III instrument components

Beam delivery in the YAG laser head

A Nd:YAG laser (1) (mini-laser module) generates the therapy beam. The beam passes through the energy attenuator (5), is then combined with the aiming beam (6) and coupled into the viewing path of the corneal microscope by means of a coupling mirror (11).

The prism head of the slit illumination system (13) is lowered by 10°. The lens system (2) can be used to adjust the focus shift.

Note for VISULAS YAG III:

In the central illumination position, a small part of the therapy beam is clipped.

Note for VISULAS YAG III Combi:

The prism head of the slit illumination system must be lowered by about 30° before the YAG laser can be fired and the YAG laser beam can arrive at the patient's eye without encountering any obstacles.



Fig. 5 LSL YAG III beam delivery

Design of the VISULAS YAG III Combi

The VISULAS YAG III Combi may consist of the following components:

- 1 VISULAS YAG III laser head
- 2 LSL YAG III Combi laser slit lamp
- 3 VISULAS 532s laser console control panel
- 4 VISULAS 532s laser console
- 5 VISULAS YAG III laser console
- 6 IT 4L instrument table (optional)
- 7 Foot switch (optional)


Fig. 6 VISULAS YAG III Combi instrument components

Beam delivery in the LSL YAG III Combi laser slit lamp

The laser beam is interlaced into the illumination beam path, so that the slit image and the laser spot can move around the fundus in synchrony.

This ensures that the area of the fundus which is to be coagulated is always illuminated, even in the peripheral fundus.

This means of combining illumination and laser beams is particularly advantageous for small illumination fields.

- 1 Patient's eye
- 2 Contact lens
- 3 Optical system for slit lamp and laser spot visualisation
- 4 Corneal microscope with physician's safety filter
- 5 Laser beam expander
- 6 Coupling lens
- 7 Optical fibre
- 8 Slit illumination system

Fig. 7 Beam delivery in the LSL YAG III Combi laser slit lamp

Controls

VISULAS YAG III control unit


The control unit (2) is the central power supply and control unit for the VISULAS YAG III and the laser slit lamp.

The control panel (5), the LASER STOP button (4) and the combined control knob/button (3) are located on the front side of the control unit. The VISULAS YAG III is menu operated via the buttons on the control panel and the combined control knob/button. The combined control knob/button is used to adjust the system parameters in the individual operating modes by turning the knob. Pressing the button selects the *energy* parameter in all modes.

The recessed grip on the top (6) is used for safe transport of the control unit.


The key switch (1) is located on the right-hand side.

The power supply and control connectors required for operation of the VISULAS YAG III are located on the rear panel.

- 1 Key switch
- 2 Control unit
- 3 Control knob/button
- 4 LASER STOP button
- 5 Control panel
- 6 Recessed grip

Fig. 8 VISULAS YAG III control unit

- 1 Foot switch connector
- 2 Door interlock connector
- 3 Serial port for servicing
- 4 Connector for combined operation with the VISULAS 532s
- 5 Slit lamp connector (control cable)
- 6 Slit lamp connector (high-voltage cable)
- 7 Power inlet connector with integrated fuses

Fig. 9 Control unit/connectors

Note

The cables and the corresponding sockets and ports are colour-coded for easy identification. The sockets and ports are also labelled.

Both connectors and sockets are marked with red dots. These must be lined up when connecting the cables.

- To disconnect the cables, grasp the connector by the sleeve and pull it directly backwards. Pulling at the sleeve releases the lock mechanism on the connector.

Caution

Only connect cables and connectors that are the same colour and size. Never force connections.

- Install the remote interlock contact. If you do not wish to install an interlock, connect the yellow coding plug provided in place of the interlock (yellow connector, 2, Fig. 9).

VISULAS 532s control panel for VISULAS YAG III Combi

The VISULAS YAG III Combi is menu operated using the buttons on the VISULAS 532s control panel (1) and the combined knob/button (6).


To operate the laser console, which is mounted under the table on the right, the control panel is removed and placed on the table (see also the VISULAS 532s user manual).

- Insert the connector (2) into the port (3).
- To remove, hold the control panel at the front with both hands. Tilt the upper edge forwards so that the guide rod (7) is released from the retaining element (4).
- Remove the control panel by pulling it towards you carefully.

The control panel's lower guide rod (8) can be rotated by up to 90°, allowing the control panel to be safely angled towards the user.

- You can vary the system parameters in the individual operating modes by turning the combined control knob/button (6).

Pressing the button selects the *laser power output* parameter in all modes.

- 1 Control panel
- 2 Connector
- 3 Connector for control panel
- 4 Sprung retaining elements
- 5 Fixed retaining elements
- 6 Combined control knob/button
- 7 Upper guide rod
- 8 Lower guide rod
- 9 Cable holder

Fig. 10 Removing the VISULAS 532s control panel

Laser slit lamp

- 1 Laser head
- 2 Focus shift control
- 3 Fastening screw
This screw fastens the binocular tube or other component to the microscope body.
- 4 Interpupillary distance scale
- 5 Binocular tube
- 6 Eyepieces
The standard equipment includes eyepieces with a magnification of 10x. An eyepiece with a reticle can be used for taking measurements and to aid focusing.
- 7 Magnification changer
- 8 Cover for halogen lamp
(To replace lamp, see page 62)
- 9 Slit length control
Variable in steps of 1 / 3 / 5 / 9 / 14 mm.
Special slit 1 x 5 mm, $\pm 45^\circ$, 90°
- 10 Slit width control (0 to 14 mm)
- 11 Filter selector
White line pointing up:
Full aperture combined with a heat-absorbing filter.
White line to the patient:
Blue filter for tonometry and fluorescence observations.
White line to the physician:
Red-free filter for contrast enhancement during fundus observation.
- 12 Prism head
- 13 VISULAS 532s laser spot size controller

Note

A tonometer can be mounted to the slit lamp.
The tonometer support is mounted on the joint on the carrier arm. The tonometer itself is attached by a magnet and fitted only when required.


Fig. 11 VISULAS YAG III laser slit lamp


Fig. 12 LSL YAG III Combi laser slit lamp

Instrument base with headrest

- 1 Fixation lamp mounting

Caution

The laser slit lamp's fixation lamp is only to be operated with the standard red blinking diode provided. Otherwise, the fixation lamp may overheat and be damaged.

- 2 Patient headrest
- 3 Height adjustment for chinrest
- 4 Patient's chinrest
- 5 Locking screw for instrument base
This screw is used to adjust the tension of the XY motion of the instrument base.
- 6 Laser connection cable
This cable is inserted into the sockets (5/6, Fig. 9) on the control unit.
- 7 Mounting holes
Provided for mounting the instrument base to the support using the appropriate screws.
- 8 Fast-action base lock
Lever tilted towards the patient: The base is locked.
Lever tilted towards the physician: The base can be moved freely.
- 9 Brightness control for slit illumination
- 10 Joystick
X/Y fine adjustment by tilting the joystick
(X/Y coarse adjustment is achieved by shifting the complete instrument base).
- 11 Hand switch
The laser can be fired using either the hand switch or the foot switch (special accessory). The control unit menu is used to switch between the hand and foot switches.

Note for VISULAS YAG III Combi:

The VISULAS YAG III Combi has an additional micro-manipulator on the hand switch (see VISULAS 532s user manual).

- 12 Index mark for mid-position
Indicates the centre point of the height adjustment range.
The height is adjusted by turning the joystick knob.


Fig. 13 Instrument base with headrest

Intended use

The VISULAS YAG III is designed for ophthalmologic use including posterior capsulotomy and peripheral iridotomy.

Note for VISULAS YAG III Combi:

The VISULAS YAG III Combi instrument system consists of the VISULAS 532s laser system, the VISULAS YAG III laser system and the LSL YAG III Combi laser slit lamp as a general applicator for the two laser systems.

For the intended use of the VISULAS 532s, the VISULAS 532s user manual must be observed.

Warning

Users in USA:

The use of the system is subject to the Intended Use as described in the following:

This device is intended for use primarily by physicians and health care workers and may only be used under the supervision of a physician. This device will not be sold to the general public.

The VISULAS YAG III is used in ophthalmic applications, which includes posterior capsulotomy and peripheral iridotomy.

The VISULAS 532s in the VISULAS YAG III Combi is intended for ophthalmologic use including photocoagulation of ocular tissues.

The laser energy is delivered via either transpupillary delivery or intraocular endoprobe delivery.

Warning

Before firing the laser, make sure that the aiming beam is clearly visible in the treatment area!

Possible complications

Retinal bleeding and/or burns may occur particularly with excessive laser power.

Contact lenses

Always use the correct contact lens for the intended application.

Warning

Always hold the contact lens at right angles to the laser beam. Incorrect positioning of the contact lens causes distortion of the laser focus. This can lead to unwanted side effects.

Never fire the laser if you are unable to see the aiming beam in the target area.

Principle of Nd:YAG laser treatment

The VISULAS YAG III is an ophthalmic therapy laser system for non-contact incision of ocular tissue.

The radiation is generated by a Q-switched Nd:YAG laser, which emits pulses of typically 4 ns duration and a maximum energy of approximately 10 mJ at a wavelength of 1,064 nm. At the focus of the laser beam, a plasma is created which is at a high pressure and temperature, but is spatially confined. This phenomenon is termed optical breakdown or photodisruption. The plasma expands, thereby cooling rapidly, such that there is no thermal effect on the tissue.

However, the mechanical stress due to the shock wave accompanying the expansion of the plasma results in an incisional effect.


Fig. 14 Principle of photodisruption using the Nd:YAG laser

Contraindications

The VISULAS YAG III should not be used in cases of retinal detachment or cystoid macular oedema, or with patients with intraocular glass lenses. Glass lenses can be destroyed by the effect of the optical breakdown.

The following complications have been reported after Nd:YAG laser treatment:

<i>Complication</i>	<i>Incidence</i>
Cystoid macular oedema	Estimated: < 1 %
Retinal detachment *	Total population: 1.4 % After capsulotomy: 2.8 %
Endophthalmitis	Estimated: < 1 %
Raised intraocular pressure	Estimated: < 1 %
Rubeosis iridis (diabetic)	Estimated: < 1 %

The frequency and severity of complications depends on the energy applied. The lower the energy used, the lower the degree of postoperative inflammation and increase in intraocular pressure.

Warning

The VISULAS YAG III may only be operated by physicians who have adequate knowledge of the medical application of the instrument, the effect of the laser beam on tissue and possible side effects.

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Klin. Monatsbl. Augenheilkd. 1994; 205:181-186

© 1994 F. Enke Verlag Stuttgart

Installation and transport of the VISULAS YAG III

The VISULAS YAG III will be installed and operated for the first time by a Carl Zeiss Meditec authorised specialist. If you wish to use the VISULAS YAG III as a mobile system, you will require the optional transport box.


Caution

The VISULAS YAG III should be transported in the transport box provided only.

Before transporting the instrument, the following steps should be taken:

- Disconnect all cable connections.
- Disconnect the binocular tube from the slit lamp.
- Unscrew the screws fixing the slit lamp base to the instrument table.
- Remove the key from the key switch on the control unit.
- Adjust the slit lamp to the lowest height using the joystick.
- Unscrew the fixation lamp.
- Put all components into the cardboard boxes as shown in Fig. 15. Reassemble the equipment in reverse order. The connections are shown in Fig. 9.

- 1** Transport case lid
- 2** Transport case handle
- 3** Transport case (720 mm x 550 mm x 370 mm)
- 4** Transport case, open
- 5** Foot switch *
- 6** Laser slit lamp with headrest
- 7** Control unit
- 8** Transport castors
- 9** Fixation lamp
- 10** Accessories, packed
- 11** Laser safety goggles *
- 12** Binocular tube
- 13** Contact lenses *
- 14** Arm rests

* optional equipment, not included as standard

Fig. 15 VISULAS YAG III in transport case

Installation of the VISULAS YAG III Combi

The instrument will be installed by a Carl Zeiss Meditec authorised specialist.


Caution

Optical fibres and application systems from other manufacturers are not to be used with this laser system.

- Avoid impacts and blows during transport of the instrument.
- Mount the system in a suitable location close to the physician and the slit lamp.

To install the laser console on an IT 4L instrument table, the table mounting (1, Fig. 16) must first be fitted. Four holes are provided for this purpose under the tabletop (on the right hand side, seen from the physician's perspective).

- Screw the table mounting's mounting plate firmly to the table top using the screws provided.
- Now affix the laser slit lamp to the table.
- Slide the VISULAS YAG III from the front onto the upper plate of the table mounting.

1 Table mounting

Fig. 16 IT 4L with table mounting

Caution

To avoid damage to the instrument, remove the key before carrying out this procedure.

- Finally, holding it vertically, slide the VISULAS 532s, from the front, onto the lower plate of the table mounting.
- Place the foot switch in a suitable position under the table and insert the connector into the red socket on the bottom of the VISULAS 532s (see VISULAS 532s user manual).


Note:

The cables and the corresponding sockets and ports are colour-coded for easy identification. The sockets and ports are also labelled.

Connectors and sockets are both marked with red dots. These must be lined up when connecting the cables.

- To disconnect the cable, grasp the connector by the sleeve and pull it directly backwards. Pulling at the sleeve releases the lock mechanism on the connector.

1 VISULAS YAG III laser console

2 VISULAS 532s laser console

Fig. 17 IT 4L with laser consoles

**Caution**

Only connect cables and sockets that are the same colour and size. Never force connections.

- Install the remote interlock contact. If you do not wish to install an interlock, connect the yellow coding plug provided in place of the interlock (yellow connector, see VISULAS 532s user manual).
- The power cable for the IT 4L should be connected to a fixed power socket.
- Connect the VISULAS YAG III to one of the sockets on the instrument table (see separate IT 4L user manual). The VISULAS 532s must be plugged into a separate fixed earthed socket.
- Remove the protective cap from the optical fibre of the applicator.

**Caution**

Never touch the end of the optical fibre with your fingers. Protect the end of the optical fibre from contamination. Because of the laser's high power output, any dirt may lead to damage to the end surface of the fibre.

- Insert the optical fibre carefully into the socket for the optical fibre on the VISULAS 532s console (see VISULAS 532s user manual) and screw in firmly.


**Caution**

Screw in by hand only. Do not use tools!

- Fold out the control panel foot, unwind as much of the control panel cable as required and connect it to the front of the instrument (see also VISULAS 532s user manual).
- Adjust the foot to the required angle and place the control panel in the location required.

**Note:**

An additional component (e.g. an optical splitter) can be mounted between the binocular tube and the microscope body. This assembly is locked by the locking screw (1).

- 1 Locking screw
- 2 Binocular tube
- 3 Eyepieces

Fig. 18 Binocular tube and eyepieces

Moving the IT 4L instrument table with the instrument

To avoid damage when moving the complete system, proceed as follows:

- Carefully lower the table to the lowest position.

Caution

Ensure that there are no objects in the path of the tabletop when moving the table downwards.

- The table can now be lifted slightly by the narrow side of the tabletop, on the slit lamp side, and moved on the castors.
- Use both hands to lift the table.

Laser treatment using the VISULAS YAG III / VISULAS YAG III Combi

Adjusting slit lamp and eyepieces

Before treating a patient with laser radiation, you should be fully familiar with the use of the slit lamp. You should therefore carefully read the user manual, and in particular the safety notes.

- Ensure that the eyepieces are inserted fully and that the eyecups are pulled out (pushed in for spectacle wearers).
- Rotate the diopter setting ring on both eyepieces fully in the "+" direction (counter clockwise).
- Place a sheet of paper (business card or similar) in the object plane of the slit lamp.
- Centre the slit projector and binocular tube to the slit lamp base.
- Select the highest magnification setting on the slit lamp.
- Open the slit.
- Look through the eyepieces. Focus on the surface of the paper using the joystick and lock the slit lamp base.
- Close the illuminating slit to a narrow line.
- Select the lowest magnification setting on the slit lamp.
- Looking through the eyepieces, turn the diopter setting ring on both eyepieces in the "-" direction (clockwise) until the surface of the sheet of paper is sharply in focus.
- Select the other magnification settings on the slit lamp.
- The image should remain in focus for all magnifications. If it does not, the procedure should be repeated.
- Make a note of the eyepiece settings for future use. For future sessions, the eyepieces must simply be adjusted to this value.

Note:

If several physicians use this instrument, it is advisable to draw up a table with the individual refractive powers and keep it in a handy location near the instrument.

Start-up and laser treatment

Control of the VISULAS YAG III is menu-driven using the control panel. Functions and menus are described in detail below.

VISULAS YAG III: Description of menus

Note for VISULAS YAG III Combi:

The menus for the VISULAS YAG III Combi are described on page 48.

Menu navigation

The VISULAS YAG III is operated via a menu system with four menu screens via the VISULAS YAG III control unit (laser console).


Fig. 19 VISULAS YAG III menu navigation

Power-up mode

When the system is switched on using the key switch (1, Fig. 8) on the VISULAS YAG III control unit, the starting screen appears on the control panel. The instrument is now in **Power-up mode** (Fig. 20).

Fig. 20 Power-up mode

The system automatically performs a self-test. The starting screen includes a progress bar showing the progress of the system test. The following actions are carried out and the following tests and initialisation routines run during the self-test:

- The system is switched on and initialised
- The system monitor is started
- The watchdog circuit is tested and then activated
- The energy power down is tested
- Internal laser shots are fired and energy values determined for each pulse mode

When the system test has been successfully completed, the program switches into **Treatment STANDBY mode** (Fig. 22).

You can switch directly into **Diagnosis mode** during the system test by touching the **ECO** button.

Diagnosis mode

Diagnosis mode (Fig. 21) enables the user to use the laser slit lamp as a diagnostic slit lamp.

The VISULAS YAG III control unit is used merely as a power source for the slit lamp.

The laser remains switched off and the hand switch and foot switch are deactivated.

Fig. 21 Diagnosis mode

Pressing the **ARROW** button, switches the program to **Power-up mode** and then automatically to **Treatment STANDBY mode**.

Treatment STANDBY mode

After the system has been switched on and the internal system test has been completed, the system enters this mode automatically unless the instrument is switched to **Diagnosis mode** by the user.


Fig. 22 Treatment STANDBY mode



All treatment-related laser parameters are displayed on the control panel and can be changed in this mode.

The following laser parameters can be set in this mode: energy, pulse mode and aiming beam intensity.

Each time the system is restarted, the energy and the aiming beam intensity values most recently selected will be used.

- To set individual laser parameters, touch the appropriate parameter button and set the value required by turning the control knob.
- Select the desired pulse mode (single, double or triple pulse) by touching the appropriate button. The reverse-coloured button (blue on a white background) is activated.
- Pressing the control knob takes you directly to the entry field for the laser energy. Rotate the knob to set the energy required.

The buttons and displays in **Treatment mode** have the following functions:

	<p>Displays the selected focus shift mode ANT - anterior focus shift POST - posterior focus shift 0 - Zero focus shift See <i>Focus shift and defocusing</i> on page 56.</p>
	<p>Laser warning lamp This field lights up in yellow when the instrument is in treatment READY mode and the aiming beam is switched on. This warning lamp lights up in red when the laser is fired.</p>
<p>Selection and display of pulse mode (burst) (here: single pulse mode activated)</p>	
	<p>Single pulse</p>
	<p>Double pulse</p>
	<p>Triple pulse</p>
<p>Selection and display of energy levels</p>	
	<p>Energy levels When this field is activated, the energy can be set using the control knob. There are 22 energy levels. In single-pulse mode, the lowest level corresponds to about 0.3 mJ, the highest to about 10 mJ. The energy display indicates the average energy of the last five shots, calculated for the currently selected energy level. The output energy of the Nd:YAG laser may vary slightly. You should therefore check the displayed energy value regularly.</p>
<p>Energy and pulse counters</p>	
	<p>Display field Displays the sum of the emitted pulses and the cumulative energy.</p>
	<p>Reset button Touching this button resets the sum of emitted pulses and the cumulative energy to zero.</p>
<p>Information mode</p>	
	<p>Switch to information mode Touching this button switches to Information mode. See page 44.</p>

Aiming beam	
	Switch aiming beam on. In Treatment STANDBY mode, the aiming beam is switched on and off by touching this button. If the aiming beam is on, this button appears in reverse colour.
	Aiming beam intensity After touching this button, the intensity of the aiming beam can be adjusted by turning the control knob.
STANDBY/READY	
	Display and switch between STANDBY/READY Touching this button switches the system into Treatment READY mode after approximately three seconds. This transition time is visually indicated by a blinking button. Touching this button again returns the system to Treatment STANDBY mode. The key illustrated is the STANDBY key.

 **Note**
If the system has been in **Treatment STANDBY mode** for more than 30 minutes, on switching to **Treatment READY mode** the system will release a series of internal laser shots in order to update the energy values for each pulse mode.

Information mode

Touching the **INFO** button in **Treatment STANDBY mode** switches the system into **Information mode** (Fig. 23).


Fig. 23 Information mode

This mode allows you to view system information. You can also set the firing mode (hand switch or foot switch), the system time (date and time), the user interface language and the control panel LC display contrast.

The following system information is displayed:

- ☐ Company name/instrument name
- ☐ Software version
- ☐ Operating time counter
- ☐ System time (date and time)
- ☐ Information on any faults which have occurred

The following options can be set in **Information mode**:

	<p>Firing mode</p> <ul style="list-style-type: none"> Touch the appropriate button to select either hand switch or foot switch firing. <p>The foot switch is not included as standard, but is available as an option.</p>
	<ul style="list-style-type: none"> Touching this button switches the system into Diagnosis mode (see page 41)
	<p>Language</p> <p>The user interface is available in the following languages: German, English, French and Spanish.</p> <ul style="list-style-type: none"> Touch the appropriate button to select the language required.
	<p>Contrast</p> <p>Adjustment of the LC display contrast to the ambient brightness</p> <ul style="list-style-type: none"> Touch this button and adjust display contrast by turning the control knob.
	<p>System time</p> <ul style="list-style-type: none"> Touching this button takes you to the system time field and allows you to set the system time (Fig. 24).


Fig. 24 Setting the system time

Setting the system time


- Activate each numerical field in turn by touching it and set the current date and time by turning the control knob.
- To save the settings and return to **Information mode**, touch the  button.
- Touching the  button also returns the system to **Information mode**, but **without** saving any previous changes.

Touching the **ARROW** button, switches the system from **Information mode** into **Treatment STANDBY mode**. The settings made are retained when the instrument is switched off.

Treatment READY mode

- Touching the **STANDBY/READY** button whilst in **Treatment STANDBY mode** switches the system, after about three seconds, into **Treatment READY mode** (Fig. 25).

1 Laser warning lamp

Fig. 25 Treatment READY mode

In **Treatment READY mode**, the aiming beam and the laser warning lamp are always on.

Parameters are set in the same way as for **Treatment STANDBY mode** (see page 42).

When the firing switch is activated, the laser beam is fired according to the parameters set.

Warning

Never fire the laser if the aiming beam is not clearly visible in the target area.

- Touching the **STANDBY/READY** button switches the system back into **Treatment STANDBY mode**.

Note

If, in **Treatment READY mode**, more than five minutes elapses without the laser being fired, the system returns to **Treatment STANDBY mode**.

VISULAS YAG III Combi: Description of menus

Note for VISULAS YAG III:

The menus for the VISULAS YAG III are described on page 39.

Menu navigation

The VISULAS YAG III is operated via a menu system with eight menu screens using the VISULAS 532s control panel.


Fig. 26 VISULAS YAG III Combi menu navigation

**Note**

Correct operation of the VISULAS 532s is essential for safe use. You should therefore ensure that you are thoroughly familiar with the VISULAS 532s user manual before using the VISULAS 532s.

To operate the VISULAS YAG III in combination with the VISULAS 532s, the two laser consoles must be connected by a special cable.

As soon as this cable is connected, the buttons on the YAG III laser console control panel will be deactivated and combination mode will be selected.

VISULAS 532s power-up mode

- ☐ When the instrument is switched on using the key switch on the VISULAS 532s console (see VISULAS 532s user manual) the start-up screen is displayed on the control panel (Fig. 27).
- ☐ The system automatically performs a self-test.
- ☐ The user interface includes a progress bar, which shows the progress of the system.
- ☐ After the self-test has been successfully completed, the system switches to the **Main menu** (see VISULAS 532s user manual).
- ☐ You can switch directly into **Diagnosis mode** during the system test.
- ☐ To do so, touch the **ECO** button.

Fig. 27 Power-up mode

Diagnosis mode

This mode enables the user to use the LSL YAG III Combi laser slit lamp as a diagnostic slit lamp. The VISULAS 532s is used merely as a power source for the slit lamp.

- ☐ The laser module and temperature regulator remain switched off, meaning that the ventilators are operated at minimum speed (minimal noise production).
- ☐ The shutter remains closed in diagnosis mode.
- ☐ Touching the **ARROW** button switches the system into **System test** mode

Fig. 28 Diagnosis mode

VISULAS YAG III power-up mode

- ❑ When the YAG system is switched on using the key switch (1, Fig. 8), the system automatically performs a self-test (Fig. 29). When the self-test has been completed, the system switches into combination mode and message Fig. 30 is displayed on the YAG console. The YAG laser can now be operated using the VISULAS 532s control panel only.
- ❑ The combination message continues to be displayed on the YAG console as long as the system remains in combination mode and there are no system errors (Fig. 30).
- ❑ Should an error occur on the VISULAS YAG III whilst in combination mode, the VISULAS 532s logs off automatically. The relevant error message is displayed on the YAG console display. The YAG laser cannot be operated or logged back on until the error has been resolved and the system restarted.
- ❑ If an error occurs during the system self-test, an error message is displayed. If this error cannot be resolved using the table on page 58, please inform our service department.

YAG III Treatment mode

- In the main menu, select the **YAG** laser as applicator.


Fig. 31 Main menu after selecting the YAG laser as applicator

Fig. 29 System test



Fig. 30 Combination mode activated

- Select the **YAG firing mode** (hand/foot).

STANDBY

- Touching the **ARROW** button switches the system from the **Main menu** into **Treatment STANDBY mode**.


Fig. 32 Treatment STANDBY mode

Selected focus shift mode display:
 ANT - anterior focus shift
 POST - posterior focus shift
 0 - Zero focus shift
 See *Focus shift and defocusing* on page 56.


Laser warning lamp:
 This field lights up in yellow when the instrument is in treatment READY mode and the aiming beam is switched on. This warning lamp lights up in red when the laser is fired.

In treatment mode, the therapy beam can be applied via the LSL YAG III Combi laser slit lamp.

All treatment-related laser parameters are displayed on the control panel and can be changed in this mode.

- To set individual laser parameters, touch the appropriate parameter button and select the value required by turning the control knob.
- Select the pulse mode required.
- Pressing the control knob takes you directly to the entry field for the laser energy of the therapy beam. Rotate the knob to set the energy required.
- The aiming beam can be turned on and off using the **AIMING BEAM** button. When the aiming beam is on, this button is highlighted.
- After touching the **BRIGHTNESS** button, the intensity of the aiming beam can be adjusted by rotating the control knob.
- A pulse counter is also displayed. This can be reset to zero by touching the **RESET** button.
- Touching the **ARROW** button returns you to the **Main menu**.
- Touching the button for switching between treatment modes **STANDBY/READY** switches the system, after about three seconds, into **Treatment READY mode**. This transition time is visually indicated by a blinking button.

READY


Fig. 33 Treatment READY mode

- ☐ In **Treatment READY mode**, the aiming beam and laser warning lamp are always on.
- ☐ Parameters are set in the same way as for **Treatment STANDBY mode**.
- ☐ When the firing switch is activated, the laser beam is fired according to the parameters set.

- ☐ Touching the **STANDBY/READY** button switches the system back into **Treatment STANDBY mode**.

- ☐ Touching the **ARROW** button returns you to the **Main menu**.

**Note**

If, in **Treatment READY mode**, more than five minutes elapses without the laser being fired, the system returns to **Treatment STANDBY mode** and emits a short beep.

Using the aiming beam system with the VISULAS YAG III / VISULAS YAG III Combi

The VISULAS YAG III uses a four-point aiming beam system. If the image is free from distortion, the following aiming beam spot patterns can be seen, depending on the position of the focus:


Fig. 34 Focusing the aiming beam

When working with energies of up to 1.5 mJ, you can focus directly on the target plane. The four spots then merge into a single spot. When working with higher energies, you must use anterior or posterior defocusing, depending on the individual situation (→ *Focus shift and defocusing*, page 56).

Where there is astigmatic distortion, the aiming beam spot pattern will take on a different appearance, usually exhibiting rhombic distortion. In this case, it is not always possible to merge the four aiming beam spots when focusing on the target plane. From this it should be clear, that in some circumstances it may no longer be possible to achieve the desired effect at low energy levels, as the therapy beam will be similarly distorted. It may be necessary to use a higher energy level in order to achieve the desired effect. Return to a normal energy level when a square aiming beam pattern is seen and it is possible to merge the four aiming beam spots.

Note

Always work with the lowest energy possible.

Intraocular lenses made of silicone are considerably more sensitive than intraocular lenses made of PMMA. For this reason, you should be particularly careful when working with silicone lenses.

The aiming beam spot pattern may also be distorted if the contact lens is held at too much of an angle. You should therefore avoid holding the contact lens at too much of an angle for the purpose of minimising reflections.

Warning

Never fire the laser if you cannot clearly see the aiming beam spot in the target area.


Fig. 35 Distorted aiming beam


Fig. 36 Focus shift control

Focus shift and defocusing

The optimal mechanical effect of the optical breakdown occurs at a specific distance from the focus of the therapy beam. For this reason, the aiming beam and therapy beam of the VISULAS YAG III are not parfocal, but are shifted (focus shift). Depending on the application, you can select one of three settings:

☐ Posterior:

Labelled as **post**, switch to the right.

The therapy beam is focused 150 μm posterior to the aiming beam focus. This setting should only be used for applications in which the area to be treated is posterior to the tissue to be protected (see also Fig. 14).

☐ Anterior:

Labelled as **ant**, switch to the left.

The therapy beam is focused 150 μm anterior to the aiming beam focus.

This setting should only be used for applications in which the area to be treated is anterior to the tissue to be protected.

☐ Zero:

Labelled as **0**, switch in mid-position.

The therapy beam and aiming beam are focused in the same plane (parfocal). This setting should only be used for applications in which there is no tissue which needs to be protected either directly anterior or posterior to the treatment region.

The laser is intended for use with the focus in these three positions only.

Note

Always turn the focus shift control until it clicks into position (an audible click should be heard).

The focus shift of the VISULAS YAG III is designed for pulse energies of up to about 1.5 mJ.

☐ If you wish to work with higher energy levels (e.g. for fibrotic secondary cataracts), you must increase the focus shift distance by defocusing:

- In switch position **ant.**, you must defocus in the anterior direction.
- In switch position **post.**, you must defocus in the posterior direction.

Note

The degree of defocusing required depends on factors such as the application and the specific conditions of the eye to be treated. It is therefore not possible to give a general recommendation.

Troubleshooting tables

Faults are indicated as system information on the control unit display. The user is also alerted by short beeps. Remedy the fault according to the troubleshooting tables below.

If a fault occurs which cannot be remedied by using the following tables, please contact a Carl Zeiss authorised specialist, citing the error message displayed.

Label the instrument as being out of order and remove the key from the key switch.

Error messages on the VISULAS YAG III control unit

<i>Message</i>	<i>Remedy</i>
LASER STOP is activated; to continue, please unlock.	<ul style="list-style-type: none"> • Unlock the LASER STOP button.
Please release actuator!	<ul style="list-style-type: none"> • Release foot or hand switch.
Door interlock is active, please close door.	<ul style="list-style-type: none"> • Shut the door or check that the shorting plug is correctly plugged into the yellow socket.
Temperature deviation, please wait!	<ul style="list-style-type: none"> • The internal monitoring systems have detected an excess temperature. Please wait a moment. • If this fault occurs frequently, ensure that the instrument's ventilation slots are not covered and that the ambient temperature falls below 35 °C.
Slit lamp not connected!	<ul style="list-style-type: none"> • Check that the slit lamp is properly connected.
No foot switch connected!	<ul style="list-style-type: none"> • Connect the foot switch or operate the instrument using the hand switch. Select the required firing mode in the <i>Information mode</i>.
In the STANDBY and READY modes, demo is displayed beside the info button. The laser cannot be fired.	<ul style="list-style-type: none"> • Switch the instrument off and then on again. While the progress bar is active during the system test (see <i>Power-up mode</i>, page 40), touch the ZEISS logo on the screen, then press the control knob/button.

An error has occurred on the VISULAS YAG

Faults during operation of the VISULAS YAG III Combi

If, in YAG mode, an error occurs on the VISULAS YAG III, an error message will be displayed on the control panel.

At the same time, an error message with a specific error code will be displayed on the YAG console.

Remedy the fault according to the table on page 57, then switch the VISULAS 532s off and then on again.

Faults on the VISULAS 532s are indicated as system information on the control unit display. The user is also alerted by short beeps.

Note: Follow the instructions given.

Warning message: Remedy the fault according to the troubleshooting tables below.

If a fault occurs, which cannot be remedied by using the following tables, please contact a Carl Zeiss Meditec authorised specialist, citing the error number displayed.

Label the instrument as being out of order and remove the key from the key switch.

Error messages and information on the VISULAS 532s control unit

<i>Message</i>	<i>Remedy</i>
Door interlock is active, please close door.	Shut the door or check that the shorting plug is correctly plugged into the yellow socket.
LASER STOP is activated; to continue, please unlock.	Unlock the LASER STOP button.
No physician safety filter or slit lamp connected!	Check that the slit lamp and the OPMI (ENDO physician's safety filter) are correctly connected and properly inserted.
No fiber connected!	Check that the optical cable is properly inserted and screwed in tightly.
No foot switch connected!	Check that the foot switch is correctly connected and properly inserted.
Temperature deviation, please wait!	If this fault occurs frequently, ensure that the instrument's ventilation slots are not covered and that the ambient temperature falls below 35 °C.
Please swivel lamp housing.	Swivel the lamp casing out of position. Check that all the connectors on the rear of both consoles are inserted properly.
Pulse interval matched.	If the laser parameters 'pulse length' or 'power' are changed, the system may carry out an adjustment of the 'pulse interval' parameter. To avoid overloading the system, the interval between laser pulses is increased.
More than 20 % power deviation.	The internal monitoring systems have identified that the actual laser power output deviates from the target power output selected on the control panel by more than 20 %. Should this error message occur repeatedly, please inform our service department.
The foot switch is activated.	Release the foot switch.

<i>Message</i>	<i>Remedy</i>
Monitoring disabled, for service purpose only!	You are in a special service mode. All error messages are deactivated. Please contact our service department. Patients must not be treated in this mode!
Applicator LINK is selected, are you sure?	Safety question, to check that the correct applicator has been selected. Confirm by touching the <i>tick</i> , or abort by touching the <i>cross</i> .
Applicator LIO is selected, are you sure?	Safety question, to check that the correct applicator has been selected. Confirm by touching the <i>tick</i> , or abort by touching the <i>cross</i> .
Error 15	A button has been pressed for more than 5 s. Release the button and continue with the ARROW button.
In the STANDBY and READY modes, during YAG operation, demo is displayed beside the aiming beam button. The laser cannot be fired.	Switch the VISULAS YAG III off and then on again. While the progress bar is active during the YAG III system test, touch the ZEISS logo on the screen, then press the control knob/button.
In the STANDBY and READY modes, during 532s operation, demo is displayed beside the aiming beam button. The laser cannot be fired.	Switch to STANDBY mode, touch the demo button for more than 5 seconds.

Laser slit lamp faults

Fault	Possible cause	Remedy
No function at all	<ul style="list-style-type: none"> - Control unit power cable not connected. - Control unit not switched on. - Connecting cable not connected to control unit. - Fuse defective. 	<ul style="list-style-type: none"> - Connect the power cable. - Switch on the control unit. - Insert connector (5/6, Fig. 9). - Replace fuse. (→ page 61)
No slit illumination	<ul style="list-style-type: none"> - Slit is closed. - Halogen lamp defective. 	<ul style="list-style-type: none"> - Adjust slit width by turning control knob (10, Fig. 11). - Replace the halogen lamp (→ page 62).
Poor slit illumination	<ul style="list-style-type: none"> - Lamp voltage setting too low. - Halogen lamp not properly inserted. - Incorrect slit length setting. - Incorrect slit width setting. - Filter selector not correctly set. 	<ul style="list-style-type: none"> - Adjust halogen lamp to the required brightness by turning control knob (9, Fig. 13). - Insert halogen lamp properly (→ page 62). - Adjust slit length by turning control knob (9, Fig. 11). - Adjust slit width by turning control knob (10, Fig. 11). - Check position of filter selector (white line should be up).
Difficulties when looking through corneal microscope	<ul style="list-style-type: none"> - Corneal microscope magnification not properly adjusted. - Interpupillary distance of binocular tube not properly adjusted. - Eyepieces not properly adjusted. 	<ul style="list-style-type: none"> - Turn magnification changer knob (7, Fig. 11) to select the desired magnification. - Adjust interpupillary distance of binocular tube. - Adjust eyepieces.
Movement of instrument base too stiff	<ul style="list-style-type: none"> - Instrument base locked with fast-action base lock. - Base locked with locking screw. 	<ul style="list-style-type: none"> - Unlock fast-action base lock (8, Fig. 13). - Release locking screw (5, Fig. 13).

Replacing the fuse on the VISULAS YAG III control unit


The control unit fuses are integrated into the combined power inlet connector/fuse assembly. This assembly is located on the rear side of the control unit (→ 7, Fig. 9).

- Turn the instrument off at the power switch.
- Disconnect the power cable.
- Open flap (1, Fig. 37) using a suitable tool.
- Remove the defective fuses and insert new ones.

Note

Use fuses with the correct rating only.
(→ *Technical Specifications*, page 70)

1 Flap

Fig. 37 Replacing the fuse


Replacing the halogen lamp in the laser slit lamp

Warning

Before replacing the lamp, switch the system off at the key switch and disconnect the power cable. Allow the lamp to cool down or wear protective gloves.

- Replace the halogen lamp (3) as shown in Fig. 38.
- To unlock the lamp holder, turn the lever (5).
- Do not touch the bulb with bare fingers.
- When inserting the bulb, ensure that the lamp holder (2) is correctly positioned in relation to the locating pins (4).

- | | |
|---|--------------------------|
| 1 | Halogen lamp cover |
| 2 | Lamp holder |
| 3 | 12 V / 30 W halogen lamp |
| 4 | Locating pins |
| 5 | Lever |

Fig. 38 Replacing the halogen lamp

Care of the instrument

Warning

Disconnect the power cable before cleaning.

Cleaning optical components

The multilayer T* (T star) antireflection coating of the optical components (e.g. eyepieces and objectives) results in excellent image quality. Even slight contamination of the optics or a fingerprint will impair image quality. In order to protect internal optical elements from dust, never store the instrument with objective, binocular tube or eyepieces removed. After use, cover the instrument to protect it from dust. Always store objectives, eyepieces and accessories that are not being used in dust-free cases.

The external surfaces of optical components (eyepieces, objectives), can be cleaned as required:

- Dust on the optical surfaces can be removed using a squeeze blower or a clean, grease free brush.
- Thorough cleaning can be performed quickly and easily using a moist, anti-static cleaning cloth.
Please follow the instructions on the cleaning cloth packaging.

This will enable you to maintain the high image quality of your product.

Cleaning of painted surfaces

All painted surfaces can be cleaned with a damp cloth.

Do not use aggressive or abrasive cleaning agents.

Use spray or wipe disinfectants for cleaning and disinfecting the instrument case, foot switch, chin rest and forehead strap.

The display should be cleaned with a lightly moistened cloth only.

Ensure that no moisture penetrates the instrument during cleaning and disinfection.

Paper pads should be replaced after each patient.

Disposal

This instrument contains electronic components. At the end of its lifetime, the instrument should be disposed of in accordance with the relevant regulations.

Disposal of the product within the EU

EU directives and national regulations currently in force at the time of marketing prohibit the disposal of the product specified on the delivery note in domestic waste or by municipal waste disposal companies.

For more information on disposal of the product, please contact your local dealer, the manufacturer or their legal successors. Please also take note of the manufacturer's current online information.

If the product or its components are resold, the seller has the duty to notify the buyer that the product must be disposed of in accordance with currently valid national regulations.

Safety checks

To ensure proper functioning of the instrument, it should be subjected to regular safety inspections.

These safety inspections should be performed annually by the manufacturer or a manufacturer authorised specialist and the results recorded in the instrument logbook.

These inspections should also be performed when the instrument has been out of use for more than one year.

The inspection should include the following procedures:

<i>Measure</i>	<i>Procedure</i>
Electrical safety	Protective earth resistance with power cable: $R_{PE} < 0.2 \text{ Ohm}$. Measurements should be carried out on all conductive parts which could carry hazardous voltages in the event of a fault.
	Earth leakage current: Earth leakage current is to be measured while the instrument is in operation using a standard meter (Measurement instructions must be observed!). $\leq 0.5 \text{ mA}$ in acc. with IEC 60601-1 or (in Germany only) $\leq 1 \text{ mA}$ stand-by unit leakage current in acc. with DIN VDE 0751
	Visual inspection
Safety inspection	Foot switch, hand switch
	Safety switch-off, if actual energy deviates from nominal energy
	Remote interlock contact
	Threshold values for pulse modes 1, 2, 3
	Optical breakdown in air
	Function test
Acceptance protocol	Check adjustment of the therapy laser to the slit lamp.
	Check adjustment of the aiming beam to the therapy laser beam.
	Check internal energy meter.
Cleaning of the slit lamp lens	
Preparation of the instrument for use and handing it over to the user.	

Note for VISULAS YAG III Combi:

Safety inspection	Visual check of the physician's safety filter
	Safety switch-off, if the actual energy deviates from the nominal energy
Acceptance protocol	Check the power output of the aiming beam and therapy beam downstream of the applicator

Calibration of energy measurement

The VISULAS YAG III is calibrated such that the energy setting on the control panel coincides with the actual energy delivered to the area treated. The calibration of the energy measurement system must be checked at least once every 12 months.

The US Food and Drug Administration (FDA) demands that manufacturers of Class III and IV medical laser instruments provide US customers with energy measurement system calibration procedures.

Warning

Energy measurement system calibration is a service procedure to be performed only by trained and authorised service technicians. Because of the risk of exposure to laser radiation, this procedure is dangerous. Improper calibration can lead to serious injury.

The control unit and the laser console may only be opened by Carl Zeiss Meditec service technicians or by other persons who have been expressly authorised to do so in writing.

The manufacturer accepts no liability for damage or injury caused by unauthorised personnel tampering with the instrument. Such actions will also render any warranty claims invalid.

These instructions are written for service technicians who have been officially trained by the manufacturer.

Possession of service literature (service manual, installation instructions or similar) and/or the required tools does not entitle unauthorised persons to calibrate or otherwise service this laser device.

Requirements

Before carrying out calibration of the instrument, the optics of the laser system must be:

- ☐ perfectly adjusted and
- ☐ absolutely clean.

Alignment and cleaning of the optical system may only be carried out by authorised service technicians.

Calibration procedure

Tools and aids required for calibration of the laser system:

- ☐ Normal tools
- ☐ Energy meter for the measurement of ns pulses (calorimeter or pyrometer) with a measuring range of 0.1 to 100 mJ
The energy meter must be properly adjusted and calibrated. Before use, wait for the energy meter to attain room temperature (approx. 30 min).
- ☐ Laser safety goggles

Proceed as follows:

- Switch the focus shift control to the **post.** position.
- Remove the laser head cover.
- To do so, unscrew the two fastening screws on the top of the cover and the fastening screw of the focus shift knob.
- Remove this knob and carefully remove the cover by lifting vertically.

Warning

High voltage components are now accessible. Do not touch any components on the ignition board.

- Affix the energy meter in front of the slit lamp.
- Switch on the laser using the key switch (*laser will enter treatment **STANDBY** mode*).

Warning

All persons in the laser hazard area while the laser is on must wear laser safety goggles.

- Turn off the slit illumination.
- Set the following laser parameters:

ENERGY	maximum level
PULSE MODE	single pulse
AIMING BEAM	10
- Direct the aiming beam onto the active surface of the energy meter. Adjust it so that the square formed by the four aiming spots covers 30 to 70 % of the active surface. The energy meter must not be positioned at a distance of more than 90 mm from the objective of the slit lamp.
- Lock the base of the slit lamp.
- Set the AIMING BEAM brightness to 1 and switch to treatment **READY** mode.
- Fire five laser shots. Record the values measured by the external meter and calculate the mean of the five values.
- Compare the ENERGY value displayed on the control unit with the calculated mean. The ENERGY value displayed must not deviate from the calculated mean by more than +/- 10 %.

- Adjust the ENERGY value displayed on the control unit until it lies within the tolerance range. To do so, use the potentiometer on the square PCB at the top of the laser head (towards the binocular tube):
Turn the potentiometer:
clockwise: to increase the ENERGY
counter clockwise: to reduce the ENERGY
- Fire five laser shots and calculate the new mean. Readjust the potentiometer if necessary.
- Repeat this process until the displayed ENERGY is within the tolerance range.
- Reattach the cover and remount the focus shift control.

Ordering data, accessories and spare parts

Overview of components

	Catalogue No.
VISULAS YAG III with control unit and laser head	000000-1272-175
Laser slit lamp YAG III	000000-1272-071
Laser slit lamp YAG III Combi	000000-1172-069
Straight binocular tube (parallel tube) f = 140 mm	319770-9011-000
Optional: convergence tube f = 140 mm	319770-9001-000
Eyepiece 10x	319770-9110-000
Optional: 12.5x eyepiece	319770-9112-000

Contact lenses

Abraham capsulotomy contact lens	306877-9001-000
Abraham iridotomy contact lens	306877-9011-000

Laser safety goggles

Safety goggles YAG in acc. with DIN EN 207 (IR 1064 L5); 190 to 400 nm; 1064 nm	000000-0450-052
Laser safety goggles Argon / 532 in acc. with DIN EN 207, 400 to 532 nm	000000-0450-051

Fuses

2x T 3.15 A / E, 5 x 20 mm in acc. with IEC 60127	000000-0302-915
---	-----------------

Foot switch

Foot switch, complete	000000-1149-630
-----------------------	-----------------

Transport case

Transport case, complete	000000-1285-536
--------------------------	-----------------

Laser training eye

Laser training eye with holder (with 10 spare capsules, membranes, fundi and irises)	000000-1122-486
--	-----------------

Miscellaneous

Halogen lamp, 12 V 30 W for slit lamp	000000-0120-704
Face mask	301496-0000-000
Table mounting for VISULAS YAG III Combi on IT 4L	000000-1268-416

Carl Zeiss Meditec and their authorised representatives will be happy to provide you with information on further accessories.

- ☐ DigiCam adapter
- ☐ Video documentation equipment
- ☐ Applanation tonometer

VISULAS YAG III laser system

Laser class	IV (CFR 21, Part 1040, Sec. 1040.10) 4 (DIN EN 60825-1:2003)		
Wavelength of therapy beam	1064 nm		
Mode	Supergaussian		
Pulse length	typ. 4 ns		
Pulse mode	Energy (typical)	Max. pulse repetition rate	Burst frequency
1 (single pulse)	10 mJ	2.5 Hz (5 pulse / 2 s)	-
2 (double pulse)	23 mJ	1 Hz (1 pulse / s)	33 kHz
3 (triple pulse)	35 mJ	0.5 Hz (1 pulse / 2 s)	33 kHz
Energy attenuation	22 levels: 2, 4, 6, 8, 10, 12, 14, 16, 20, 24, 28, 32, 36, 40, 42, 48, 56, 60, 64, 70, 80 and 100 % transmission		
Beam diameter at the focus	10 µm in air (1 / e ²)		
Angle of exit aperture (divergence)	16° (round angle)		
Aiming beam	Wave length: 660 to 680 nm Power output: max. 150 µW 4-point aiming beam system for focusing		
NOHD	2 m		
Rated voltage	100 V to 240 V ± 10 %		
Rated frequency	50 to 60 Hz		
Rated current	Max 1.4 A to 0.7 A; transient: max. 2.8 A (< 1 s)		
Fuses	2x T 3.15 A / E, 5 x 20 mm in acc. with IEC 60127		
Electrical Protection Class	SK I		
Protection type	IP 20		
Instrument type	B (in acc. with IEC 60601-1)		
Earth conductor	The instrument should only be connected to properly earthed power outlets.		
Control unit dimensions	H x W x D = 135 mm x 210 mm x 330 mm		
Control unit weight	4 kg		

LSL YAG III laser slit lamp

Motion range of instrument base	Lateral: 110 mm Horizontal: 90 mm Vertical: 30 mm
Slit lamp	12 V 30 W halogen lamp, adjustable. Slit length variable in steps of: 1 / 3 / 5 / 9 / 14 mm; Slit width continuously adjustable: 0 to 14 mm Slit image rotation: 0°, ±45°, 90°
Corneal microscope	Magnification with magnification changer: 5x, 8x, 12x, 20x, 32x with 10x eyepieces and f = 140 mm tube. Straight binocular tube, f = 140 mm with PD adjustment from 55 to 78 mm. Optional: convergence tube Optional: 12.5x eyepieces
Dimensions (incl. laser head)	H x W x D = 625 mm x 300 mm x 450 mm
Weight (incl. laser head, tube and eyepieces)	11 kg
Electrical power supply	Electrical power is supplied by the VISULAS YAG III laser console.
Protection type	IP 20

LSL YAG III Combi laser slit lamp

Motion range of instrument base	Lateral: 110 mm Horizontal: 90 mm Vertical: 30 mm
Slit lamp	12 V 30 W halogen lamp, adjustable. Slit length variable in steps of: 1 / 3 / 5 / 9 / 14 mm; Slit width continuously adjustable: 0 to 14 mm Slit image rotation: 0°, ±45°, 90°
Corneal microscope	Magnification with magnification changer: 5x, 8x, 12x, 20x, 32x with 10x eyepieces and f = 140 mm tube. Straight binocular tube, f = 140 mm with PD adjustment from 55 to 78 mm. Optional: convergence tube Optional: 12.5x eyepieces
Dimensions (incl. laser head)	H x W x D = 623 mm x 350 mm x 400 mm
Laser treatment spot size	Continuously adjustable from 50 to 1000 µm (without contact lens) parfocal, larger spot sizes depending on contact lens used
Laser beam delivery	Interlaced with slit illumination system
Weight (incl. laser head, tube and eyepieces)	12.0 kg
Electrical power supply	Electrical power is supplied by the VISULAS YAG III laser console.
Protection type	IP 20
Accessories	Tonometer, co-observation tube, video documentation equipment etc. from the range of accessories for the SL 120 and SL 130 slit lamps

Ambient conditions for overall system

Ambient conditions for overall system	Ambient temperature: 10 to 40 °C Relative humidity: 0 to 90 % (non-condensing) Air pressure: 700 to 1060 hPa
Transport and storage conditions	Ambient temperature: -25 to 55 °C Relative humidity: 0 to 90 % Air pressure: 700 to 1060 hPa

The VISULAS YAG III complies with the EC Medical Device Directive 93/42/EEC and the national implementation of this directive in the form of the German Medical Products Act (MPG).

Device class under MPG: II b

UMDNS No.: 16-947 (VISULAS YAG III)

UMDNS No.: 12-281 (laser slit lamp)

This declaration is rendered invalid if changes are made to the product without the manufacturer's authorisation.

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CARL ZEISS MEDITEC AG